---

中文版 | English

VITLAB®

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VITLAB®

## micropipette -8/-12

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---

## Contents

1.Safety Instructions	4
2.Purpose	5
3.Limitations of use	5
4.Operating limitations	5
5.Operating exclusions	5
6.Operating and Control Elements	6
7.Pipetting	7
8.Checking the Volume	9
9.Accuracy Table	10
10.Adjustment	11
11.Autoclaving	12
12.UV sterilization	12
13.Servicing and Cleaning	13
14.Disassembling the Manifold	15
15.Ordering Information • Accessories	19
16.Spare Parts	19
17.Troubleshooting	20
18.Repairs • Calibration Service	21
19.Warranty Information	22
20.Disposal	22

---

---

# 1. Safety Instructions

## Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe all specifications provided by reagent manufacturers.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 5). If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visible damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 20), and contact the manufacturer if needed.

---

## 2. Purpose

Air-displacement pipette for pipetting aqueous solutions of medium density and low to medium viscosity.

## 3. Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C to 40 °C. Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 50 kPa
- Viscosity: 260 mPa · s

## 4. Operating Limitations


Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than  $\pm 1$  °C

## 5. Operating Exclusions

The user has to ensure the compatibility of the instrument with the intended application. This instrument cannot be used:

- for liquids incompatible with polypropylene and FKM
- for liquids incompatible with polycarbonate (inspection window)
- for liquids that corrode polyvinylidene fluoride and silicone
- for liquids that corrode polyphenylene sulfide (on 50 µL, 100 µL, 200 µL and 300 µL instruments)
- for liquids of a very high vapor pressure

# 6. Operating and Control Elements

(Fig. shows VITLAB® micropipette 8 channels)

\* Volume display  
The figures in the display are read from top to bottom, the dash represents the decimal point.

**Note:**  
Optimum analysis results can only be obtained with quality tips. We recommend BRAND® pipette tips.

# 7. Pipetting

## 1.Fitting the tip

Use the correct tips according to the volume range or the color code.  
Ensure that the tips are securely seated.  
Pipette tips are disposable items!


## 2.Volume setting

Select the desired volume by rotating the volumes setting wheel.  
Avoid twisting and abrupt rotating motions during this adjustment.

## 3.Align the manifold

The manifold can turn freely in both directions.

## 4.Aspirate sample

- a)Press pipetting button to the first stop.
- b)Hold the pipette vertically and immerse the tip 2 to 3 mm into the liquid.

Volume range	Immersion depth in mm	Waiting time in s
> 0.5µl -100 µL	2 -3	1 -2
> 100 µl -300 µL	2 -4	1 -2

- c)Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tip immersed for another 1-2 s.
- d)Touch the tip against the container wall.

#### 5. Discharge sample

- a) Place the pipette tips against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.
- b) Press the pipetting button slowly to the first stop and hold it down. For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- c) The blow-out stroke empties the tips completely: Press the pipetting button down to the second stop.
- d) While doing this, wipe the pipette tips against the wall of the container.
- e) Remove the pipette tips from the container wall and let the pipetting button slide back.

#### 6. Ejecting the tips

Hold the manifold over a suitable disposal container and press the tip ejection key to the stop.

#### **Note:**

ISO 8655 prescribes rinsing the pipette tips once with the sample liquid prior to the actual pipetting process.

#### **Important!**

Don't lay the instrument horizontal when the tips are filled. Liquid may enter and contaminate the instrument.

## 8. Checking the Volume

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

### 1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 7 for procedure).

### 2. Condition the pipette

Condition the pipette before testing by using a pipette tips to aspirate and discharge the test liquid (distilled H<sub>2</sub>O) five times. After this, discard the pipette tips.

### 3. Carry out the test

a) Attach new pipette tips and prerinse one time with test liquid.

b) Aspirate liquid and pipette it into the weighing vessel.

#### Note:

Each individual channel must be tested separately.

c) Weigh the pipetted quantity with an analytical balance. (Please follow the operating manual instructions from the balance manufacturer.)

d) Calculate the volume, taking the temperature into account.

e) At least 10 pipettings and weighings in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

### Calculation for nominal volume $V_0$

$x_i$  = Weighing results

$n$  = Number of weighings

Mean value

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume

$$\bar{V} = \bar{x} \cdot Z$$

Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

$Z$  = Correction factor (e.g. 1,0029 µL/mg at 20 °C, 1013 hPa)

Accuracy\*

$$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

Coefficient of Variation\*

$$CV\% = \frac{100 s}{\bar{V}}$$

\*Calculation of accuracy (A%) and variation coefficient (CV%): A% and CV% are calculated according to the formulas for statistical control.

#### Note:

Testing instructions (SOPs) are available for download at [china.brand.com.cn](http://china.brand.com.cn)

# 9. Accuracy Table

VITLAB® micropipette -8/12, Digital adjustable

Volume range μL	Volume step μL	A* ≤ ± %	CV* ≤ %	Increment μL	Type of tip μL
0.5- 10	10	1.6	1.0	0.01	0.5-20
	5	2	2		
	1	8	6		
5 - 50	50	0.8	0.4	0.05	2-200
	25	1.4	0.8		
	5	6	3		
10 - 100	100	0.8	0.3	0.1	2-200
	50	1.4	0.6		
	10	4	2		
20 - 200	200	0.8	0.3	0.2	2-200
	100	1.4	0.6		
	20	4	1.5		
30 - 300	300	0.6	0.3	0.5	5-300
	150	1.2	0.6		
	30	3	1.5		

Final test values related to the nominal capacity (maximum volume) indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C) and with smooth operation. According to DIN EN ISO 8655.

\* A = Accuracy


\* CV = Coefficient of Variation

---

## 10. Adjustment

The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made.

1. Check the volume, determine actual value (see page 9).
2. Remove the cover: Push the hook forward, raise it slightly and then pull it back.
3. Using a paperclip or an unused pipette tip, remove the protective film (this protective film can be discarded).

4. Push the red adjustment slider completely back (right arrow), raise the volume-setting wheel (left arrow) and release the adjustment slider.

5. Set the adjustment value:  
set the volume-setting wheel to the previously determined actual value. A volume check is recommended after every adjustment.

6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Reinsert the cover.

---

## 11. Autoclaving

The pipette is completely autoclavable at 121 °C, 200 kPa absolute with a holding time of at least 15 minutes, according to DIN EN 285.

1. Eject the pipette tips.
2. Autoclave the complete pipette without any further disassembling.
3. Allow the pipette to completely cool and dry.

### Note:

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

### Attention!

Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between).

If the pipette is autoclaved frequently, the pistons should be oiled with silicone oil in order to preserve smooth movement. Please use only the recommended silicone oil, see accessories page 19. If necessary after sterilization, tighten the connection between the hand grip and the manifold.

## 12. UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

---

## 13. Servicing and Cleaning


The three main components of the manifold can be easily separated and disassembled for servicing, cleaning or replacing parts. The procedures are illustrated on the following pages.

**Note:**

The changing of the O-rings on the individual nose cones is described in detail in the instructions enclosed with the spare part.

### 13.1 The main components of the manifold

- A Piston unit with piston support bar [ I ] and pistons inserted in this unit, which can be unscrewed individually for cleaning or replacement.
- B Nose cone assembly with nose cone support bar [ II ] and central guide rod (Z), which is attached to this, and the nose cones and seals, which can be unscrewed individually for cleaning or replacement.
- C Manifold housing, which is connected to the manifold housing cover [ III ] of the piston unit with two turnlock fasteners.

---

## 13.2 Servicing

In order to assure proper functioning, the VITLAB® micropipette -8/-12 should be serviced and cleaned at regular intervals.

What is to be inspected?

- a) Inspect nose cones, pistons and seals for damage and contamination.
- b) Test the sealing of the instrument. To do this aspirate a sample, and then hold the instrument in a vertical position for about 10 sec. If a drop forms at the tip orifices, see the troubleshooting guide, page 20.

## 13.3 Cleaning instructions

- a) Clean single nose cones, pistons and nose cone support bar/ piston support bar (these components only) with soap solution or isopropyl alcohol. Afterwards, rinse with distilled water.
- b) Let these parts dry and cool down completely. Residual moisture in the nose cones may result in a loss of accuracy.
- c) Lubricate the piston with a very thin coating of the silicone oil supplied. For the central guide rod (Z) only use the recommended fluorstatic grease! See accessories page 19.

## 13.4 Disconnecting handle from manifold

- a) Eject the pipette tips.
- b) To disconnect the manifold, pull it downward as far as possible, and only then turn it clockwise. After one rotation, it should no longer be pulled downward while it is being turned.

### Attention:


Improper handling can damage the unit!


---

## 14. Disassembling the Manifold


### 14.1 Removing of nose cones and seals


- 1) Disconnect the manifold from the handle.
- 2) Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.


- 3) Push the mounting tool (M) (supplied with the pipette) on the nose cone and unscrew it.


- 4) Push the piston unit down to the bottom. The seal will either remain inside the nose cone or will stay on the piston after the nose cone is removed. Remove the seal, inspect it and clean or replace if necessary. If required, lubricate the piston lightly with the supplied silicone oil.

(In the VITLAB® micropipette -8/-12 30-300  $\mu$ L unit, remove the additional pressure ring from the piston for cleaning.)

- 5) Push the seal on the piston with its flat side facing upward. (In the VITLAB® micropipette -8/-12 30-300  $\mu$ L unit, first replace the additional pressure ring!)


Use the mounting tool to mount and tighten the cleaned or new nose cone.


- 6) Reassemble the manifold. The manifold must be screwed into the handle counter-clockwise, until it audibly snaps into place. Check the instrument for tightness (see page 14, 'Servicing').

## 14.2 Removing of pistons for cleaning or replacing


- 1) Disconnect the manifold from the handle.
- 2) Turn both closures of the manifold housing cover 90° (e.g., using a coin) and slide off the housing.


- 3) Remove both outer Phillips screws on the nose cone assembly.

**Attention:**


Do not loosen the central guide rod (Z)!


- 4) Pull the piston and nose cone assembly apart and separate.  
Remove stroke springs.


- 5) Place the mounting tool on the piston nut and unscrew the piston nut.

6) Remove the piston nut with piston spring and pull out the piston.

7) Insert the cleaned or new piston and piston spring ( a ). Screw on the piston nut and tighten with the mounting tool. Replace the stroke springs ( b ).

8) Reassemble the manifold. The manifold must be screwed counter-clockwise into the handle so that it audibly snaps into place.

Check the instrument for tightness, mobility and accuracy (see page 14, 'Servicing').

# 15. Ordering Information and Accessories

## 15.1 Ordering Information

VITLAB® micropipette -8

Volume	Cat. No.
0.5 - 10 µL	C0708900
5 - 50 µL	C0708906
10 - 100 µL	C0708908
20 - 200 µL	C0708910
30 - 300 µL	C0708912

VITLAB® 12 micropipette -12

Volume	Cat. No.
0.5 - 10 µL	C0708920
5 - 50 µL	C0708926
10 - 100 µL	C0708928
20 - 200 µL	C0708930
30 - 300 µL	C0708932

## 15.2 More Accessories

Bench-top rack for 6 VITLAB® micropipettes or 3 VITLAB® micropipettes -8/-12  
Cat. No. VT003.1672002

Silicone oil for VITLAB® micropipette  
5 mL/10 mL Cat. No. 703677

Fluorstatic grease Cat. No. 703678

Shelf/rack mount for 1 VITLAB® micropipette bzw. VITLAB® micropipettes -8/-12  
Cat. No. VT003.1672000

Reagent reservoir, non-sterile, PP,  
Volume 60 mL, VE 10 Cat. No. 703459

# 16. Spare Parts

Parts will differ slightly depending on nominal volume of instrument.  
(Fig. shows spare parts for VITLAB® micropipette -8/-12, 10-100 µL).

Capacity	A	B*	C	D
0.5 - 10 µL	705659	705677	703380	703340
5 - 50 µL	705666	705634	705618	703343
10 - 100 µL	705662	705635	705618	703344
20 - 200 µL	705663	705636	705618	703345
30 - 300 µL	705664	705638	705618	703346

\* incl. seal, O-ring and mounting tool.  
VITLAB® micropipette -8/-12 30-300 µL  
with additional pressurering.

## 17.Troubleshooting

Problem	Possible cause	Corrective actio
Tip dripping (instrument leaks)	- Unsuitable tip - Tip not seated tightly	- Only use high-quality tips - Press tip on firmly
The instrument does not aspirate or aspirates too little; the discharged volume is too low	- Seal contaminated - The seal or cone is damaged - The piston is contaminated or damaged	- Clean seal - Replace seal or shaft - Clean or replace piston
Aspiration is too slow	- Shaft clogged	- Clean shaft
Discharged volume is too large	- Pipetting button pressed too far into the blow out position before sample uptake	- Operate properly. See "Pipetting" page 7
Piston is difficult to move	- The piston is contaminated or needs oil	- Clean and oil piston
Tip dripping (instrument leaks)		

---

## 18. Repairs - Calibration Service

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!

### 18.1 Return for repair

- a) Clean and decontaminate the instrument carefully.
- b) Complete the 'Declaration on Absence of Health Hazards' (ask your supplier or manufacturer for the form).
- c) Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used.

The return transport of the instrument is at risk and cost of the sender

### 18.2 Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

The detailed testing instruction can be downloaded on [china.brand.com.cn](http://china.brand.com.cn)

We also offer a paid calibration service for your equipment. Please contact the manufacturer or distributor to obtain a BRAND CNAS laboratory calibration request form. After completing the calibration request form and the 'Declaration on Absence of Health Hazards', simply send the equipment to be calibrated and the completed form to us. The calibration laboratory will return the calibration report along with the equipment within fifteen working days, following the standard procedures.

---

## 19. Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

## 20. Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice.

We will not be held responsible for printing or typographical errors.